

DermAware


ISTRUZIONI PER L'USO

band with velcro fastener




1 WEARING

press for
turning ON
white led
turns ON




press for
turning OFF
white led
turns OFF




2 ON / OFF


Level 1
Vibration + Blue LED



Level 2
Vibration + Green LED




Level 3
Vibration + Red LED



15 seconds
cycle

3 FUNCTIONALITY




Arduino code upload through programming board
connected via USB



Re-programming through Wi-Fi connection
with Arduino

4 PROGRAMMING / RE-PROGRAMMING



Charge through
Micro USB port

Remove the bottom screws
to replace the battery

5 CHARGE/ REPLACE BATTERY

Please carefully read the above mentioned instructions
before trying the DermAware device.

While wearing the DermAware smart bracelet, it is
normal to feel an uncomfortable experience: it is part of
the process hearing an unpleasant sound and feeling a
stressful vibration.

Please do not wet the device and store it in a dry
environment.

In case of malfunctions and different experiences than
the usual vibration and sounds, please remove the
bracelet and contact the assistance.

6 WARNINGS