


Palpatine AR: How to use it



INSTRUCTIONS:

STEP 1: Download the Palpatine APP and print out the Marker in the Marker section;

STEP 2: Place the adhesive Marker in the middle of the Bra, following the guidelines already present on the bra;



STEP 3: Wear your Palpatine bra and position your tablet on the support!
Place your self in front of the tablet.



STEP 4: Open the App on your tablet and follow the instructions to effectuate the self-palpation process;

STEP 5: Remember to repeat the process once a month!

STEP 6: Spread the word about Palpatine AR to increase the awareness about the importance of Self Palpation.

